



In-Flight Laboratory Analysis

September 2012

David Baumann, HHC Element Manager, JSC

Gail Perusek, In-Flight Lab Analysis Project Manager, GRC

Emily Nelson, PhD, Project Scientist

Michael Krihak, PhD, ExMC Lab Analysis Project Lead, ARC

Dan Brown, Systems Engineering Lead, ZIN Technologies



Overview



FY11 Highlights

- Space Environment
- Human Research Program
- Risks & Functional Requirements



Space Environment



- Zero gravity
- Radiation
- Limited mass, volume and power
- Limited resources - water, air, food
- Communication lags or blackouts (isolation)




Clay Anderson centrifuges blood samples for a nutrition project during Increment 15



ISS030-E-257690 (26 April 2012) --- European Space Agency astronaut Andre Kuipers, Expedition 30 flight engineer, prepares for IMMUNE venous blood sample draws in the Columbus laboratory of the International Space Station.



HRP Organization & Elements



- The Program is divided into 6 major elements, which
 - Provide the Program's knowledge and capabilities to conduct research, addressing the human health and performance risks
 - Advance the readiness levels of technology and countermeasures to the point of transfer to the customer programs and organizations
- The National Space Biomedical Research Institute (NSBRI) is a partner with the HRP in developing a successful research program.



Human Research Program



Program Goals

- Perform research necessary to understand and reduce spaceflight human health and performance risks in support of exploration
- Enable development of human spaceflight medical and human performance standards
- Develop and validate technologies that serve to reduce medical risks associated with human spaceflight



Human Health and Countermeasures (Research)



Research Objectives

- Developing knowledge, capabilities, countermeasures, & technologies to mitigate highest risks to crew health & performance & enable human space exploration.
- Developing technologies to reduce medical and environmental risks, and to reduce human systems resource requirements.
- Requirements align to current gaps and tasks as defined by the IRP that are research and technology development related.
- As new knowledge is gained, the required approach to research and development may change.
 - Gaps and tasks are periodically updated, subject to change.
 - *Hence, requirements are not comprehensive and allow for flexibility as new research needs are identified and arise.*



ExMC Risk and Gap (Med Ops)



Exploration Medical Capability

- Ensure astronaut health and safety due to injury or illness on extended (>30 days) human exploration missions.
- Provide biomedical diagnostics capability to facilitate the recognition and treatment of several medical conditions.
 - Address ExMC gap 4.05: Lack of minimally invasive inflight laboratory capabilities with limited consumables required for diagnosing
- Provide analysis capability of biological fluids (i.e. blood, urine, saliva, sweat) in any habitable location



Levels of Care



Level of Care	Mission	Capability
I	*LEO < 8 days	Space Motion Sickness, Basic Life Support, First Aid, Private Audio, Anaphylaxis Response
II	LEO < 30 days	Clinical Diagnostics, Ambulatory Care, Private Video, Private Telemedicine
III	Beyond LEO <30 days	Limited Advanced Life Support
IV	Lunar > 30 days	Medical Imaging, Sustainable Advanced Life Support, Limited Surgical, Dental Care
V	Mars Expedition	Autonomous Advanced Life Support and Ambulatory Care, Basic Surgical Care

*LEO = Low Earth Orbit (e.g. Space Shuttle orbit)



In-Flight Laboratory Analysis



Specifications

- Minimize the equipment's mass, volume, consumables, reagents and power.
- Ease of operation; minimal training.
- Ideally, should have FDA approval, or have gone through the rigors of FDA approval type validation.
- As mission duration lengthens, an analyzer's capability should be readily expanded through software, reagents, dipsticks and/or microfluidic cartridges.
- Short start-up time
- 3- to 5-year shelf life



Functional Requirements



Technology Measurement Capability

Element	Analytes (Examples)	Biological Sample	Technological Approaches
Radiation	Damage markers Genetic, epigenetic traits of individual susceptibility	Cell cultures Blood Blood cells	Lab-on-a-chip (LOC) Microscopy Immunofluorescence (microscopy, flow cytometry)
Behavioral Health and Performance	Fatigue, stress markers Circadian rhythm markers Depression markers	Blood	LOC biomarker detection
Human Health and Countermeasures	Bone resorption, metabolism markers Ca homeostasis hormones Natural Ca isotope composition Phosphorus, potassium, magnesium Blood cell surface markers, cytokines Oxidative stress markers Bone metabolism and hormonal regulation markers Biomarkers: Physiologic, Metabolic, Oxidative Stress, Cardiovascular	Blood Urine Sweat Cell models - <i>in vitro</i> Cell culture	Point-of-care devices LOC diagnostics Flow cytometry Mass spectroscopy Colorimetric and fluorimetric microplate analysis Biosensors



Classes of Research Measurements

Analyte Class	Examples
Ions (Na, Cl, etc)	Na, Cl, K,...
Blood Gases	pH, pO ₂ , pCO ₂ , BUN,...
Small Molecules	Glucose, lactate,...
Amino acids	3-methylhistidine, GABA,...
Proteins	Il-1, leptin, transferrin, troponin,...
Peptides	BNP, helical peptide P, insulin...
Enzymes	ALT, AST, CK-MB,...
Fatty Acids	Triglycerides,...
Minerals	Fe, Zn, Se, Cu, Mg, P,...
Vitamins	Retinol, b-carotene, folic acid..
Steroids	Cortisol, estradiol, DHEA,...
Lipids	Cholesterol, LDL, HDL,...
Metabolites	Bilirubin, creatinine,...
Cell Type	Leukocyte, WBC, hematocrit,...
Cell Markers	P-selectin, CD4,...



Functional Requirements



Operational Requirements

MEASUREMENTS

Basic Metabolic Panel	Blood Gases Panel	Hematology	Cardiac Panel	Liver/ Renal Panel	Urinalysis
Glucose Calcium Sodium Potassium CO ₂ , Total Chloride BUN Creatinine	PaO ₂ PaCO ₂ SaO ₂ HCO ₃ pH	WBC Count RBC Count HCT Hgb Neutrophils Abs. Neutrophils Count Lymphocytes Monocytes Eosinophils PLT	Troponin I CK-MB	Total Bilirubin Direct Bilirubin ALP AST ALT	Specific Gravity pH Leukocytes Nitrites Proteins Glucose Ketones Urobilirubin Bilirubin Blood Urate



Instrument Requirements



- OPERATION
 - Start-Up Time - 2 minutes
 - Completion of Analysis - give indication
 - Self-Calibration
 - User - single caregiver or patient (ease of use)
- SHELF-LIFE
 - Storage
 - Cartridge storage up to 3 years at room temperature
- POWER
 - Battery - operate up to 144 hours
 - Switchable power source - external or internal sources



One-Year Evaluation



Expectations

- Perform a multiplexed analysis of whole human blood (demonstrated using serum is acceptable) of the following four analytes:
 - (25 OH) Vitamin-D
 - N-terminal telopeptide (NTx)
 - IFN- γ (interferon gamma)
 - TNF- α (tumor necrosis factor alpha)
- Demonstrate assay performance against a gold standard
- Demonstrate multiplexed assay (4-plex)
- Demonstration at JSC



One-Year Evaluation



Expectations (continued)

Measurement Ranges

- NTx:
 - Adult Male: 5.4-24.2 nM BCE;
 - Premenopausal, Adult Female: 6.2-19.0 nM BCE
- TNFalpha: 0-22 pg/mL
- IFNgamma: 0-5.0 pg/mL
- Vitamin D: 20-150 ng/mL

(BCE – bone collagen equivalents)



Summary



- One-year study objectives align with HRP requirements
- HRP requirements include measurement panels for research and medical operations
 - These measurement panels are distinctly different.
- Instrument requirements are defined
 - Power, volume and mass not quite a critical limitation as for medical operations (deep space exploration missions)
- One-year evaluation goals will lead HHC towards in-flight laboratory analysis capability